

## How to send Primescam / Omnicam cases to Image Gallery

### Add Image Gallery as a lab in Sirona Connect:

Log onto Sirona Connect (via the internet)

<https://customer.connectcasecenter.com/>

Go to **"My Account"**

Click on **"My Favorite Laboratories"**

Click on **"Search Labs"** (towards bottom of screen)

Fill in Country and Zip Code

Click on **"Start Search"**

Scroll down until you see **"Image Gallery"** and click the plus sign to Add them to your favorite labs

When you log onto Sirona Connect to send a case, Image Gallery will be listed.

### Sending cases to Image Gallery

#### ADMINISTRATION PHASE:

**Restoration type:** select **"IMPLANT"** (or other restoration)  
**Design Mode:** choose **"SCREW RETAINED"** OR **"CEMENT RETAINED"**  
**Scan body type:** for **"MANUFACTURER"**, select manufacturer of implant  
for **"IMPLANT PLATFORM"**, select implant  
for **"SCAN BODY"**, select **"ATLANTIS IO FLO"**

#### ACQUISITION:

Follow scanning prompts as usual

#### CONNECT PORTAL:

**"CHECK RESTORATION DATA"** – Confirm the information is correct  
**"ENTER ORDER DATA"** – Fill out Patient Information. Uncheck Anonymize if you want the patient name connected to the order. If you check anonymize an order number will be sent instead of patient name.  
**"ADDITIONAL FILES"** – Attach any photos or other files here  
**"RECIPIENTS"** select **"IMAGE GALLERY"**  
**"ADDITIONAL ORDER INFORMATION"**: add any extra notes to the lab here  
**"RETURN DATE / TIME"**  
**"SUBMIT CART"**: a new screen populates and shows you which cases have been "commissioned".