

## Porcelain Fused to Metal (PFM) - Anterior

*Noritake*

<u>Ingredients (Porcelain)</u>	<u>Percent</u>
Tin Oxide (as glass)	< 26.0
Silicon Dioxide (as glass)	< 51.6
Aluminum Oxide (as glass)	< 11.5
Calcium Oxide (as glass)	< 0.6
Magnesium Oxide (as glass)	< 0.5
Potassium Oxide	< 7.0
Sodium Oxide (as glass)	< 7.4
Lithium Oxide (as glass)	< 0.3
Pigments	< 1.1

### Ingredients (Metal)

**IDENTALLOY<sup>®</sup>**  
CERTIFICATE

HIGH NOBLE ALLOY CONTAINS AT LEAST 40% GOLD AND 60% (TOTAL) GOLD, PALLADIUM AND PLATINUM.

The manufacturer certifies that the dental casting alloy provided to the laboratory with this certificate is a High Noble alloy.  
The laboratory certifies that High Noble alloy was used to fabricate this prosthesis.

Alloy	Manufacturer	Composition
Encore	DENTSPLY Ceramco	Au 48.0% Pd 40.0% Pt 0.0%

DENTIST CERTIFICATE (attach to patient record)  
[www.identalloy.org](http://www.identalloy.org)

**IDENTALLOY<sup>®</sup>**  
CERTIFICATE

HIGH NOBLE ALLOY  
24844

<b>ALLOY</b> Encore
<b>MANUFACTURER</b> DENTSPLY Ceramco
<b>COMPOSITION</b> Au 48.0% Pd 40.0% Pt 0.0%

LAB COPY