

How to send CARESTREAM (CS3600) cases to Image Gallery

After acquiring scan data, you have two options to send us files, “**CSCConnect**” and “**Send To**”

1. **CSCConnect:** Login with email and password, choose Image Gallery under “Pick Lab” or input our lab’s email address if this is the first case: chris@imagegallerydublin.com
2. **Send To:** Gives you the option of saving the file to the desktop where you can then attach the files in an email. Send to: chris@imagegallerydublin.com

If any issues arise, please contact CARESTREAM tech support @ 1-866-724-6317.